



WHITE PAPER

Opentrons Electronic Pipettes

Technical specs, test methods, and data validation

Written by
Opentrons



Table of Contents

SECTION 1

GEN2 Pipettes

Introduction	4
Single-Channel Technical Data	6
8-Channel Technical Data	7
Testing Methods	8
Single-Channel Example Gravimetric Testing Output	9
8-Channel Example Gravimetric Testing Output	11
Single-Channel Test Data	12
8-Channel Test Data	13
Technical Drawings	14
Lifetime Testing	16

SECTION 2

GEN1 Pipettes

Introduction	18
Testing Methods	21
Test Data	23
Citations	26

SECTION 1

Single and 8-Channel GEN2 Pipettes

INTRODUCTION

Single-Channel and 8-Channel GEN2 Pipette

GEN2 Pipette Product Description

Opentrons GEN2 Pipettes are designed specifically for lab automation and optimized for use with the OT-2. The Opentrons Single-Channel GEN2 pipettes reliably conduct liquid transfers from 1 μL -1,000 μL , while the Opentrons 8-Channel GEN2 pipettes reliably conduct liquid transfers from 1 μL -300 μL . Both have tip pickup and drop +/- 1mm of a successful seal when using [Opentrons tips](#).

Opentrons GEN2 pipettes feature comparable specs to pipettes 10x as expensive. They also allow researchers to complete protocols using fewer pipettes due to a newly broadened volume range and any combination of two of these pipettes.

GEN2 SINGLE-CHANNEL AVAILABLE RANGES

- P20 GEN2 (1 μL – 20 μL)
- P300 GEN2 (20 μL – 300 μL)
- P1000 GEN2 (100 μL – 1000 μL)


GEN2 8-CHANNEL AVAILABLE RANGES

- P20 GEN2 (1 μL – 20 μL)
- P300 GEN2 (20 μL – 300 μL)



FIGURE 1

Pipette Design Cutaway Diagram with Parts Labeled



Single-Channel GEN2 Pipette

TABLE 1

Single-Channel GEN2 Pipette Specifications

MODEL	VOLUME	ACCURACY		PRECISION	
		% D	µL	% CV	µL
P20 GEN	1 µL	+/- 15%	0.15 µL	+/- 5 %	0.05 µL
	10 µL	+/- 2%	0.2 µL	+/- 1 %	0.1 µL
	20 µL	+/- 1.5%	0.3 µL	+/- 0.8 %	0.16 µL
P300 GEN2	20 µL	+/- 4%	0.8 µL	+/- 2.5%	0.5 µL
	150 µL	+/- 1%	1.5 µL	+/- 0.4%	0.6 µL
	300 µL	+/- 0.6%	1.8 µL	+/- 0.3%	0.9 µL
P1000 GEN2	100 µL	+/- 2%	2.0 µL	+/- 1%	1 µL
	500 µL	+/- 1%	5.0 µL	+/- 0.2%	1 µL
	1000 µL	+/- 0.7%	7.0 µL	+/- 0.15%	1.5 µL

8-Channel GEN2 Pipette

TABLE 2

8-Channel GEN2 Pipette Specifications

MODEL	VOLUME	ACCURACY		PRECISION	
		% D	μL	% CV	μL
P20 GEN2	1 μL	+/- 20%	0.2 μL	+/- 10%	0.1 μL
	10 μL	+/- 3%	0.3 μL	+/- 2%	0.2 μL
	20 μL	+/- 2.2%	0.44 μL	+/- 1.5%	0.3 μL
P300 GEN2	20 μL	+/- 10%	2.0 μL	+/- 4%	0.8 μL
	150 μL	+/- 2.5%	3.75 μL	+/- 0.8%	1.2 μL
	300 μL	+/- 1.5%	4.5 μL	+/- 0.5%	1.5 μL


TESTING METHODS

Single and 8-Channel GEN2 Pipette

Gravimetric testing for the Opentrons Pipette GEN2 is holistic in that we test the entire system rather than just the individual pipette. Specifically, the Pipette GEN2 is tested with the OT-2 Robot, Opentrons Tips, and Opentrons API exactly as a customer would use it.

Further, Opentrons adheres to a slightly stricter version of ISO8655 for pipette gravimetric testing. The key difference is that Opentrons does not include a pre-wet test for each individual aspirate of the pipette. A graphical breakdown of the process is seen in Figure 2.

FIGURE 2
Gravimetric Calibration Test Process Diagram



EXAMPLE GRAVIMETRIC TESTING OUTPUT

Single-Channel GEN2 Pipette

Example Output from Gravimetric Testing

The standard output of this process is a series of measurements taken at the pipette's minimum and maximum volumes. The data gathered is used to generate systematic and random errors for each pipette.

Pipette Serial No: P20SV202019101717			Time/Date: 10/31/2019 10:48:00		
Min Vol CV Spec:	5	Min Vol %D Spec:	15	Scale Serial No:	NB-A-552550
Mid Vol CV Spec:	1	Mid Vol %D Spec:	2	Temperature:	22.3 C
Max Vol CV Spec:	0.8	Max Vol %D Spec:	1.5	Humidity:	58% RH

TABLE 3

Pipette Calibration Testing

TEST VOLUME (µL)	TRIALS (µL)									
	1	2	3	4	5	6	7	8	9	10
1	1.052	1.0693	1.006	0.9899	0.982	1.052	0.982	1.0073	1.0599	1.00
10	9.948	10.078	10.016	9.988	10.04	10.028	10.012	10.0733	10.0866	10.0386
20	20.07	20.078	20.054	19.97	19.98	19.936	20.0473	20.092	20.068	19.976

TABLE 4

Pipette Calibration Testing

TEST VOLUME (μL)	MEAN	STANDARD DEV	%CV	%D	RESULT
1	1.0201	0.0343611	3.369	2.0067	PASS
10	10.0309	0.0429114	0.428	0.3087	PASS
20	20.0271	0.0556160	0.287	0.1357	PASS

EXAMPLE GRAVIMETRIC TESTING OUTPUT

8-Channel GEN2 Pipette

Example Output from Gravimetric Testing

The standard output of this process is a series of measurements taken at the pipette's minimum and maximum volumes. The data gathered is used to generate systematic and random errors for each pipette.

Pipette Serial No: P20MV202020033004		Time/Date: 4/30/2020 17:41:00	
Min Vol CV Spec:	10	Min Vol %D Spec:	20
Max Vol CV Spec:	0.8	Max Vol %D Spec:	2.2
		Temperature:	22 C
		Humidity:	49.79

TABLE 5

Pipette Calibration Testing

TEST VOLUME (µL)	TRIALS (µL)									
	1	2	3	4	5	6	7	8	9	10
1	1.0273	1.0420	1.0340	1.0140	1.0560	1.0280	0.9920	1.0293	1.0360	1.0267
20	20.0220	19.9733	19.9733	19.9720	20.0040	19.936	19.9980	20.0000	19.9920	19.9907

TABLE 6

Pipette Calibration Results




TEST VOLUME (µL)	MEAN	STANDARD DEV	%CV	%D	RESULT
1	1.0285	0.01691	1.6439	2.8530	PASS
20	19.9909	0.01579	0.0790	-0.0454	PASS

TEST DATA

Single-Channel GEN2 Pipette

TABLE 7

GEN2 Single-Channel Pipette Volumetric Test Data





TEST DATA


8-Channel GEN2 Pipette

TABLE 8





GEN2 8-Channel Pipette Volumetric Test Data



Single-Channel GEN2 Pipette



8-Channel GEN2 Pipette



LIFETIME TESTING

Single and 8-Channel GEN2 Pipette

Opentrons performs lifetime testing on each pipette model to ensure customers will have reliable results for the lifetime of their pipette. Each pipette model is tested to 650+ hours of use: this equates to using 55,000 tips over 6 months with the robot pipetting 5 hours a day, 5 days a week.

There are two primary areas of wear Opentrons has addressed through lifetime testing. The first is the o-ring seal inside the pipette, and the second is the pipette nozzle that interacts with disposable tips. The o-ring seal is tested to 20-30 km of travel depending on the pipette model. The nozzle-to-tip interaction is tested for 100,000+ tip pickups and drops for each model. Gravimetric testing is performed at each $\frac{1}{4}$ life interval to ensure that the pipette is still in-spec.

SECTION 2

Single and 8-Channel GEN1 Pipettes

INTRODUCTION

Single and 8-Channel GEN1 Pipette

The OT-2 Electronic Pipettes are the first of their kind. They easily couple with the OT-2 robot's gantry, so the robot can move them precisely in XYZ space while actuating the pipette motors (more on this below).

FIGURE 3

Single-Channel Electronic Pipette



FIGURE 4

8-Channel Electronic Pipette




Electronic Pipette Design

Opentrons electronic pipettes take the classic micropipette design and simply replace the scientists' thumbs with electric motors. We use a robust stepper motor attached to a precision-ground linear shaft. The shaft creates a seal with an o-ring to create a piston (or volume chamber)

which accurately controls the volume of liquid aspirated and dispensed inside the disposable tips. The tip cones are made from a hard polypropylene plastic designed to pick up and eject 100,000 tips over the lifetime of the pipette.

FIGURE 5

Pipette Design Cutaway Diagram with Parts Labeled



Volume Calculation

The accuracy of our liquid transfers is dependent on the conversion of millimeters traveled by the pipette lead screw into a corresponding change in fluid within the pipette tip.

The first conversion that needs to happen is between the motor's microsteps and the linear motion of its lead screw. Stepper motors divide a full rotation into a series of "steps" that are moved to sequentially in order to turn the lead screw. Each step of the motor's rotation therefore translates into a corresponding travel distance of the lead screw.

Because the lead screw is connected directly to the piston in the volume cylinder, we can translate the mm traveled to the volumetric change at the tip. This is the same relationship as is found in standard piston-based manual pipettes. The correct mm to μl conversion was determined empirically by collecting volume data on each pipette at different volumes with gravimetric analysis as described below in the methods sections. Once enough data was collected, we applied a curve to the data to find the appropriate function for each pipette model.

After determining the microliters per mm function for a given pipette model, the 'constant volume' script was executed to test both the precision (%CV) and accuracy (%D) of the pipette, advertised in the Testing Data section of this whitepaper

Single and 8-Channel GEN1 Pipette

Gravimetric Capabilities and Principles

Gravimetry—the measurement of weight—is a well-accepted standard and recognized in the scientific literature as a good methodology for assessing liquid handling performance [1]. It is also relatively simple to carry out and only requires an analytical balance, an instrument found in most labs.

In using gravimetry we're taking advantage of the straightforward mass-to-volume ratio of water—that is, 1 mL of water weighs approximately 1 gram [2]. Therefore, water volume can be measured reliably on an analytical balance with the appropriate sensitivity.

Limits to Gravimetric Analysis

Gravimetry is limited by the precision of the analytical balance being used and the environmental conditions around it. When measuring small volumes, tiny variations in the lab environment like those caused by OT-2 robot's movements or a person walking close to the balance will lead to skewed results in normal laboratory settings. Additionally, gravimetry is affected by a wide variety of other environmental conditions, including evaporation, static electricity, vibration, temperature, relative humidity, and more [3].

It is important to note that, as liquid volumes become smaller, both vibrational and environmental effects become more pronounced. This is why, to ensure accurate measurement, we have limited our gravimetric analysis to volumes from 1 μL to 10 μL .

GRAVIMETRIC TESTING INSTRUMENTATION AND MATERIALS:

- Analytical balance with a 0.01 mg precision and USB connection (Radwag, AS82/220.R2)
- Low retention pipette tips (Eppendorf Low-Retention, 022493004)*
- OT-2 with modified deck
- Optical table by ThorLabs
- Windows 10 laptop computer with Microsoft Excel
- Analytical grade water (Corning, 46-000-CV, Lot 27017005)
- Draft protection (to ensure the scale does not provide false readings from the air moving in the testing room)

We have developed a custom rig that allows us to collect many gravimetric data points quickly while keeping measurement conditions the same across time.

The Z-Stage head design from our OT-2 robot was mounted on a standard 80/20 aluminum extrusion. The pipettes are attached in the same manner as they are to the OT-2 robot so that they could perform the same linear motion done during a typical aspirate or dispense command. A customized chamber made from acrylic was placed around the rig to limit airflow through the fixture, which could otherwise interfere heavily with testing volumes less than 100 μL . The rig was also placed on an optical table to dampen vibrations.

In order to reduce evaporation effects, the surface area of the water was reduced to a size of a Falcon cap. Since static electricity is another factor that affects accuracy in measurement, we used specialized tips that are less hydrophobic (Eppendorf Low-Retention, 022493004)* and a new tip was used for each reading.

The measurements were recorded from an analytical balance with a 0.01 mg precision (Radwag, AS82/220.R2) in a spreadsheet on a Windows 10 PC connected over USB 2.0 to the scale.

Gravimetric Experimental Method

Using the pipette being analyzed, water was aspirated and dispensed in place sequentially onto an analytical balance with a 0.01 mg precision (Radwag, AS82/220.R2). After waiting 2 seconds for the balance to level out after each aspirate and dispense, a measurement was recorded in a spreadsheet on a Windows 10 laptop computer attached to the scale over USB. For each measurement, 10 readings were taken sequentially. For a Multi-Channel Pipette, each channel was tested individually using the same methodology.

Data from the analytical balance was averaged and normalized automatically with the constant volume test script [4] developed by Opentrons test engineers to account for random noise in the balance and fluctuations due to environmental conditions such as evaporation.

Disclaimer

Given our updated validation process, we now recommend using Opentrons tips instead of Eppendorf. If you'd like to use other tips, please contact our contact our sales team to ask about protocol modification.

TEST DATA

Single and 8-Channel GEN1 Pipette

Gravimetric Data Processing

For our gravimetric tests, the test script extracted the measured weights from the analytical balance and output the values in a CSV file format that was then imported into an Excel spreadsheet. The weights were then converted directly into volumes according to the direct relationship between the density and volume of water, and averaged to find the test mean.

Precision and Accuracy Calculations

The calculated volumes from both gravimetric and photometric testing were then averaged and compared with the intended volumes. We assessed random error, or precision, using the following equation:

$$CV = \frac{\sigma}{X} * 100$$

CV is the coefficient of variance (expressed as a percentage).

σ is the standard deviation of the sample set.

X is the mean of all of the sample volumes.

We assessed systematic error, or accuracy, using the following equation:

$$d = \left[\frac{(\bar{X} - V_{test})}{V_{test}} \right] * 100$$

d is the systematic error (expressed as a percentage).

X is the mean of all of the sample volumes.

V_{test} is the specified test volume.

TABLE 9

Precision and Accuracy of OT-2 Electronic Pipettes at Target Volumes

PIPETTE	TARGET VOLUME (μL)	RANDOM (%CV)		SYSTEMATIC (%d)	
		%	μL	%	μL
P10 Single Channel	10	+/- 1	+/- 0.1	+/- 2	+/- 0.2
	5	+/- 3	+/- 0.15	+/- 5	+/- 0.25
	1	+/- 5	+/- 0.05	+/- 15	+/- 0.15
P50 Single Channel	50	+/- .4	+/- 0.2	+/- 1	+/- 0.5
	25	+/- 0.6	+/- 0.15	+/- 1.5	+/- .375
	5	+/- 5	+/- 0.25	+/- 5	+/- .25
P300 Single Channel	300	+/-0.3	+/- 0.9	+/- .6	+/- 1.8
	150	+/- 0.4	+/- 0.6	+/- 1	+/- 1.5
	30	+/- 1.5	+/- 0.45	+/- 3	+/- 0.9
P1000 Single Channel	1000	+/- .15	+/- 1.5	+/- .7	+/- 7
	500	+/- 0.2	+/- 1	+/- 1	+/- 5
	100	+/- 1	+/- 1	+/- 2	+/- 2
P10 8-Channel	10	+/- 2	+/- 0.2	+/- 3	+/- 0.3
	5	+/- 5	+/- 0.25	+/- 5	+/- 0.25
	1	+/- 10	+/- 0.1	+/- 25	+/- 0.25
P50 8-Channel	50	+/- 0.6	+/- 0.3	+/- 1.2	+/- 0.6
	25	+/- 1.5	+/- 0.375	+/- 2	+/- 0.5
	5	+/- 5	+/- 0.25	+/- 5	+/- 0.25
P300 8-Channel	300	+/- 0.5	+/- 1.5	+/- 1.5	+/- 4.5
	150	+/- 0.8	+/- 1.2	+/- 2.5	+/- 3.75
	30	+/- 2.5	+/- 0.75	+/- 5	+/- 1.5


Data Discussion

When compared with other liquid handling robot performance data freely available online [5], the data above shows that the OT-2 performs as well as robots 10x more expensive. Our novel electronic pipette designs have proven capable of the precision and accuracy required in executing delicate biological experiments in common laboratory environments.

Supplementary Data

Below you can find the raw data from the accuracy and precision tests run in the Opentrons lab. These links are continually updated with additional data as it becomes available.

- [P10 Single Channel](#)
- [P10 8-Channel](#)
- [P300 Single Channel](#)
- [P300 8-Channel](#)



CITATIONS

Single and 8-Channel GEN1 Pipette

1. Batista, E., Pinto, L., Filipe, E. and Van der Veen, A. (2007). Calibration of micropipettes: Test methods and uncertainty analysis. *Measurement*, 40(3), pp.338-342
2. US Geological Survey. (2017). The USGS Water Science School. [online] Available at: <https://water.usgs.gov/edu/density.html> [Accessed 23 Mar 2018]
3. TTE Laboratories (2015). ISO 8655 Pipette Standards. [online] Available at: <http://www.ttelaboratories.com/TTE-University/ISO-8655> [Accessed 23 Mar 2018].
4. Opentrons Labworks, Inc (2018). Opentrons Gravimetric Test Script. [online] Available at: https://github.com/Opentrons/Volumetric_Testing_Scripts [Accessed 25 Mar 2018]
5. Eppendorf AG (2008). EPMotion 5070 Technical Data. [online] Available at: https://online-shop.eppendorf.us/US-en/eshopdownload/downloadbykey/101139_Operating-Manual_186 [Accessed 26 Mar. 2018]*